

Account # _____
 Phone _____
 Bill to _____

Authorized By: _____
 Fax # _____
 Ship To _____

Purchase Order # _____
 Account Contact _____
 Email Address _____

Date _____
 Client Name _____
 Client Weight _____

STANDARD FEATURES

- Heavy duty frame and parts.
- 8" Casters
- 20" Rear Wheels
- 45° of tilt in space (22.5° anterior/22.5° posterior)
- Height Adjustable Stroller Handle
- Lifetime warranty on base and 5 year no question asked warranty on all Gunnell designed products.
- 90 Days on All Padding

K0009 PPC-001\$4,100.00
 Positioning Platform with no sides
 Positioning Platform with sides (SELECT PLATFORM SIDE OPTIONS)

FACTORY SETTINGS (Reference Drawing)

COLORN/C
 Bright Yellow Dark Maroon Gloss Black Red
 Medium Blue Silver Vein Black Wrinkle
 CUSTOM\$210.00
 Glitter Blue Sparkle Green Candy Red Purple Metallic
 Glitter Green Sparkle Orange Candy Blue Candy Purple
 Glitter Silver Black Stardust Green Apple Metallic

PLATFORM WIDTH (Dimension W)
 PW 27" Standard Platform Width.....N/C
 PPCM-001 " Custom Platform Width.....\$180.00

PLATFORM LENGTH (Specify Dimensions) (Dimensions A,B,C)
 PPCM-002 Platform Length
 Back (A) _____"
 Seat (B) _____"
 Leg (C) _____"
Overall Total (A+B+C) _____"
 PPCM-002 If the Overall Total is 60" or under.....N/C
 PPCM-002A If the Overall Total is greater than 60" > 70".....\$120.00
 PPCM-002B If the Overall Total is greater than 71" > 80".....\$240.00

PLATFORM TO FLOOR HEIGHT
 PFH 30" Platform to Floor Height.....N/C
 PPH-001 " Custom Platform to Floor Height.....\$370.00

HIP TO BACK ANGLE
 HBA-006 Hip to Back Angle 180° to 135°
 (Threaded Rod Adjustment Device).....N/C
 HBA-007 Hip to Back Angle 180° to 90°
 (Threaded Rod Adjustment Device).....N/C
 PPO-006 Hip to Back Angle 180° to 135°
 (Mechlok Adjustment Device).....\$780.00
 PPO-007 Hip to Back Angle 180° to 90°
 (Gas Cylinder Adjustment Device).....\$900.00

KNEE ANGLE
 PPKA-006 Knee Angle 180° to 135°
 (Threaded Rod Adjustment Device).....N/C
 PPKA-007 Knee Angle 180° to 90°
 (Threaded Rod Adjustment Device).....N/C
 PPO-008 Knee Angle 180° to 135°
 (Mechlok Adjustment Device).....\$780.00
 PPO-010 Knee Angle 180° to 90°
 (Gas Cylinder Adjustment Device).....\$900.00

PLATFORM SIDES (Dimension SH)
 PPO-001 8" Fixed Platform Sides
 (HEAD, FOOT, & SIX SIDES).....\$750.00
 PPCM-005 Fixed Sides Custom Height
 (MUST SPECIFY) _____" Desired Height
 (EACH) # _____ of sides x \$250.00\$ _____
 PPO-002 8" Flip Down Hinged Sides with Dual Locking Clamps
 (EACH) # _____ of sides x \$300.00\$ _____
 PPCM-005A Custom Height Flip Down Hinged Sides with Dual Locking Clamps
 (MUST SPECIFY) _____" Desired Height
 (EACH) # _____ of sides x \$350.00\$ _____
 PPO-002A 8" Adjustable Angle Positive Locking Sides
 (EACH) # _____ of sides x \$400.00\$ _____
 PPCM-005B Custom Height Adjustable Angle Positive Locking Sides
 (MUST SPECIFY) _____" Desired Height
 (EACH) # _____ of sides x \$450.00\$ _____


PLATFORM TILT IN SPACE
 PTIS Anterior(22.5)/Posterior(22.5) Tilt in Space.....N/C
 PPCM-006 Anterior(A)/Posterior(P) Tilt in Space.....\$400.00
 45° A / 0° P 40° A / 5° P 35° A / 10° P 30° A / 15° P
 25° A / 20° P 20° A / 25° P 15° A / 30° P 10° A / 35° P
 5° A / 40° P 0° A / 45° P

WHEEL OPTIONS
 PPRW 20" Rear Wheels.....N/C
 PPO-006 Self Propelling Gurney with 24" Wheel Drive
 with TNT Frame and No Tilt in Space.....\$320.00
 (SPECIFY SEAT TO FLOOR HEIGHT) 15" 16" 17" 18" 19"
 PPO-004 24" Wheels (ONLY available on heights 26"+).....N/C
 PPO-005 Flat Free Wheel Option.....\$125.00

STROLLER HANDLE OPTIONS
 PSH Height Adjustable Stroller Handle.....N/C
 PPO-003 Height & Angle Adjustable Stroller Handle.....\$245.00

CUSTOM OPTIONS
 PPO-003 Lateral Tilt.....\$2,500.00
 PPO-005 Tapered Corners.....\$250.00
 PPO-008 Angle Adjustable Gambled IV-pole.....\$400.00

POWER OPTIONS
 PPO-004 Power Actuation.....\$5500.00/each
 Legrests Tilt in Space Lateral Tilt Recline
 \$ _____ Total
 PPO-001 Battery Box or Vent Tray.....\$299.00
 (Specify Dimensions) Length _____" Width _____" Height _____"



FOOTNOTES ON PAGE 2

AFTER THE P.O. IS SUBMITTED ALL PLATFORM CHAIRS WILL BE CONFIRMED WITH A DRAWING FOR FINAL AUTHORIZATION. SIGNATURE REQUIRED.



BACK PLATFORM BASE PADDING OPTIONS

- PPBP-001 1" Poly Base Pad \$100.00
PPBP-002 1" Sunmate Base Pad (SPECIFY DENSITY) \$225.00
PPBP-003 2" Poly Base Pad \$185.00
PPBP-004 2" Sunmate Base Pad (SPECIFY DENSITY) \$350.00
PPBP-005 Vinyl Padding Cover with Zipper N/C
PPBP-006 Darlex Padding Cover with Zipper \$300.00
PPBP-007 Rubatex Padding Cover with Zipper \$200.00
PPBP-008 Simplex Seating System foam in place \$800.00
PPBP-009 Custom Base Padding Quote Only

SEAT PLATFORM BASE PADDING OPTIONS

- PPSP-001 1" Poly Base Pad \$100.00
PPSP-002 1" Sunmate Base Pad (SPECIFY DENSITY) \$225.00
PPSP-003 2" Poly Base Pad \$185.00
PPSP-004 2" Sunmate Base Pad (SPECIFY DENSITY) \$350.00
PPSP-005 Vinyl Padding Cover with Zipper N/C
PPSP-006 Darlex Padding Cover with Zipper \$300.00
PPSP-007 Rubatex Padding Cover with Zipper \$200.00
PPSP-008 Simplex Seating System foam in place \$800.00
PPSP-009 Custom Base Padding Quote Only

LEG PLATFORM BASE PADDING OPTIONS

- PPFP-001 1" Poly Base Pad \$100.00
PPFP-002 1" Sunmate Base Pad (SPECIFY DENSITY) \$225.00
PPFP-003 2" Poly Base Pad \$185.00
PPFP-004 2" Sunmate Base Pad (SPECIFY DENSITY) \$350.00
PPFP-005 Vinyl Padding Cover with Zipper N/C
PPFP-006 Darlex Padding Cover with Zipper \$300.00
PPFP-007 Rubatex Padding Cover with Zipper \$200.00
PPFP-008 Simplex Seating System foam in place \$800.00
PPFP-009 Custom Base Padding Quote Only

3SIDE PADDING OPTIONS

- PPSO-001 1/2" Poly Side Pads (EACH) # of sides x \$105.00
PPSO-002 1/2" Sunmate Side Pads (SPECIFY DENSITY & QUANTITY) (EACH) # of sides x \$150.00
PPSO-003 1" Poly Side Pads (EACH) # of sides x \$130.00
PPSO-004 1" Sunmate Side Pads (SPECIFY DENSITY & QUANTITY) (EACH) # of sides x \$235.00
PPSO-005 Vinyl Padding Cover with Zipper N/C
PPSO-006 Darlex Padding Cover with Zipper (EACH) # of sides x \$150.00
PPSO-007 Rubatex Padding Cover with Zipper (EACH) # of sides x \$100.00
PPSO-008 Custom Side Padding Quote Only

AFTER THE P.O. IS SUBMITTED ALL PLATFORM CHAIRS WILL BE CONFIRMED WITH A DRAWING FOR FINAL AUTHORIZATION. SIGNATURE REQUIRED.

SAMPLE AUTHORIZATION PAGE IS ON PAGE 3

DO NOT FILL OUT PAGE 3 IT WILL BE FAXED TO YOU WITHIN 48 HOURS OF RECEIVING THE P.O. FOR CONFIRMATION

FOOTNOTES

When filling out the measurement form keep in mind certain dimensions.

- 1: Inside Width: (IDW: B,S,and/or L)
Without Sides Inside Width = Outside Width of Top Platform
With Sides With Padding the Width will be Inside the Padding
With Sides Without Padding the Width will be Inside the ABS

- 1: Inside Length: (IDL)
Without Sides Inside Length = Outside Length of Top Platform
With Sides With Padding the Length will be Inside the Padding
With Sides Without Padding the Length will be Inside the ABS

Below are equations for figuring out the maximum allowable overall width

- 2: O.D. Width: (ODW)
Without Sides same as Inside Width
With Fixed Sides With Padding: 1.75" + Padding Thickness + Inside Width
With Fixed Sides Without Padding: 1.75" + Inside Width
With Hinged Flip Down Sides With Padding: 2.75" + Padding Thickness + Inside Width
With Hinged Flip Down Sides Without Padding: 2.75" + Inside Width
With Angle Adjustable Locking Sides With Padding: 3.75" + Padding Thickness + Inside Width
With Angle Adjustable Locking Sides Without Padding: 3.75" + Inside Width

3: Padding Thickness: Thickness of Both Sides of the Platform

4: Platform to Floor Height: Measured from the Ground to the Top of the Platform below the ABS

5: Maximum Fixed Sides is 8: head, foot, left side 1, left side 2, left side 3, right side 4, right side 5, and/or right side 6

6: Maximum Down or Angle Adjustable Sides 6: left side 1, left side 2, left side 3, right side 4, right side 5, and/or right side 6

7: Outside Platform Width between the Wheels on the Base: 27" with Handrims, 24" With Handrims turned Inside

8: Seat and Leg can be combined as one platform no hinge point between. Cross out leg platform, initial next to it, write full length in seat position, and cross out knee angle and initial next to it.

9: Threaded Rod Adjustment Requires Tools and is Static Adjustment


10: Mechlock or Gas Cylinders are Quick Release and are released by a cable assembly and require no tools.



GUNNELL, INC.

Dear Customer,

In an effort to serve you better, please verify your purchase order # _____ regarding the dimensions of a Gunnell Platform Chair. The sketch below describes your order as we understand it. Please verify the dimensions (change them if necessary), sign this form acknowledging your agreement and fax to (989) 871-4563. Please contact customer service at (800) 551-0055 if you have any further questions.



POSITIONING PLATFORM CHAIR: ADJUSTABLE HIP AND KNEE ANGLE

Finished Product Dimensions:

- IDLB= _____ inches Length of Back – Inside pad to hinge
- IDLS= _____ inches Length of Seat – Hip angle hinge to knee angle hinge
- IDLL= _____ inches Length of Leg Rest – Hinge to inside pad*
- IDW= _____ inches Width of Back – Inside of pads*
- IDW= _____ inches Width of Seat – Inside of pads*
- IDW= _____ inches Width of Leg Rest – Inside of pads*
- SH= _____ inches Height from top of ABS to top of sides*
- H= _____ inches Distance from top of platform steel to floor*
- ODW= _____ inches Maximum allowable cart width – for going thru doors

Fits through normal doorways Yes No

*NOTE: If no pads are ordered, finish product dimensions are inside dimensions of side walls.

PLATFORM SIDES

Number of Sides (not including head and foot): _____

- 1,4 = _____ inches Length of Back Rest Sides
- 2,5 = _____ inches Length of Seat Sides
- 3,6 = _____ inches Length of Leg Rest Sides

Options:

- Base pad thickness - _____ inches.
- Side pad thickness - _____ inches.
- Hinged sides: right _____, left _____ or both _____

The above dimensions are correct:
 _____ Date: _____

Signed _____

Print Name _____