

Account # _____
 Phone _____
 Bill to _____

 Purchase Order # _____
 Account Contact _____
 Email Address _____

 Authorized By _____
 Fax # _____
 Ship To _____

 Date _____
 Client Name _____

STANDARD FEATURES

- All Walkers
 - Heavy duty frame and parts
 - Non-Tipping Base made of durable steel
 - Height Adjustable
 - Ring Style Padding Allows for Comfortable Seating
 - Lifetime warranty on base
 - 5 year warranty on all Gunnell designed products
- Infant Walker
 1. Patient Height, Approx: 27" – 36"
 2. Balance Ring Height:
 - A. Minimum: 18"
 - B. Maximum: 24"
 3. Ring Size Opening: 8 ½" x 10"
 4. Base Width: 19"
 5. Base Length: 20"
 6. Weight (lbs): 15
- Youth Walker
 1. Patient Height, Approx: 34" – 54"
 2. Balance Ring Height:
 - A. Minimum: 22"
 - B. Maximum: 35"
 3. Ring Size Opening: 12" x 10 ½"
 4. Base Width: 31"
 5. Base Length: 34"
 6. Weight (lbs): 20
- Adolescent Walker
 1. Patient Height, Approx: 36" – 56"
 2. Balance Ring Height:
 - A. Minimum: 27"
 - B. Maximum: 42"
 3. Ring Diameter: 15"
 4. Base Width: 33"
 5. Base Length: 40"
 6. Weight (lbs): 30
- Adult Walker
 1. Patient Height, Approx: 48" – 68"
 2. Balance Ring Height:
 - A. Minimum: 43"
 - B. Maximum: 58"
 3. Ring Diameter: 16"
 4. Base Width: 34"
 5. Base Length: 47"
 6. Weight (lbs): 35

✓ SUPPORT WALKERS

- IW-1 Infant Walker.....\$559.00
- YW-1 Youth Walker.....\$749.00
- AD-1 Adolescent Walker.....\$982.00
- AW-1 Adult Walker.....\$1,199.00
- CW-1 Custom Walker.....Quote

✓ COLOR..... N/C



- | | | |
|---|---|---|
| <input type="checkbox"/> Pink | <input type="checkbox"/> Dark Maroon | <input type="checkbox"/> Purple Metallic |
| <input type="checkbox"/> Flat Black | <input type="checkbox"/> Sparkle Orange | <input type="checkbox"/> Green Apple Metallic |
| <input type="checkbox"/> Glitter Silver | <input type="checkbox"/> Glitter Blue | <input type="checkbox"/> Candy Blue |
| <input type="checkbox"/> Candy Purple | <input type="checkbox"/> Candy Red | <input type="checkbox"/> Sparkle Green |
| <input type="checkbox"/> Bright Yellow | <input type="checkbox"/> Blood Red | <input type="checkbox"/> Medium Blue |
| <input type="checkbox"/> ¹ Silver Vein | <input type="checkbox"/> Gloss Black | <input type="checkbox"/> Black Stardust |

✓ UPHOLSTRY COLORS

- Black ¹ Blue.....N/C
- Custom.....\$49.00

INFANT WALKER

E0140
YOUTH WALKER

E0140
ADOLESCENT WALKER

E0140
ADULT WALKER

E0140
REPLACEMENT PARTS

- YW-SB Scissor Board for Youth Walker.....\$99.00
- AD-SB Scissor Board for Adolescent Walker.....\$159.00
- AW-SB Scissor Board for Adult Walker.....\$199.00
- S-040 Infant Ring Pads.....\$99.00
- S-039 Youth Ring Pads.....\$109.00
- S-087 Adolescent Ring Pads.....\$119.00
- S-038-F Adult Ring Pad Rear.....\$59.00
- S-038-B Adult Ring Pad Front.....\$59.00
- S-042 Infant Sling.....\$79.00
- S-051 Youth Sling.....\$89.00
- S-088 Adolescent Sling.....\$99.00
- S-041 Adult Sling.....\$109.00

RUBATEX, FLEECE OR DARLEX COVERS

 We can cover any pad with Rubatex, Fleece or Darlex
Please call for quotes

NOTE: HCPCS Codes provided were determined by SADMERC and/or interpretations of the codes. Proper use of codes does not ensure coverage or payment. Please contact a Medicare professional for further assistance.

FOOTNOTES

1: Gunnell will select factory settings if no setting is made