

### Rehab Series Self Propelling Gurney Option

Account # \_\_\_\_\_  
 Phone \_\_\_\_\_  
 Bill to \_\_\_\_\_  
 \_\_\_\_\_  
 Purchase Order # \_\_\_\_\_  
 Account Contact \_\_\_\_\_  
 Email Address \_\_\_\_\_

Authorized By: \_\_\_\_\_  
 Fax # \_\_\_\_\_  
 Ship To \_\_\_\_\_  
 \_\_\_\_\_  
 Date \_\_\_\_\_  
 Client Name \_\_\_\_\_  
 Client Weight \_\_\_\_\_

#### STANDARD FEATURES

- Heavy duty frame and parts
- Tapered Corners (head)
- Disc Brakes with Client Operated Variable Speed Control Lever
- Front and rear anti-tips
- 2" thick Sunmate Pressure Relief foam with Darlex Cover
- 22" Platform Width
- 70" Platform Length
- Standard Platform to Floor Height with 24" wheels is 22"
- Height and Angle Adjustable Stroller Handle
- Lifetime warranty on base and 5 year no question asked warranty on all Gunnell designed products
- 90 Day Warranty on All Padding

K0009  SPG-001 .....\$4,100.00

#### FACTORY SETTINGS (Reference Drawing)

√ **COLOR** .....N/C

- Bright Yellow    Dark Maroon    Gloss Black    Red  
 Medium Blue    Silver Vein    Black Wrinkle

**CUSTOM** .....\$210.00

- Glitter Blue    Sparkle Green    Candy Red    Purple Metallic  
 Glitter Green    Sparkle Orange    Candy Blue    Candy Purple  
 Glitter Silver    Black Stardust    Green Apple Metallic

√ **PLATFORM WIDTH (Dimension IDW)**

- √ <sup>1,2</sup>SPW 22" Standard Platform Width.....N/C  
 <sup>1,2</sup>SPCW-001 \_\_\_\_\_" Custom Platform Width.....\$120.00

√ **PLATFORM LENGTH (Dimension ODL)**

- SPL Standard Platform Length 70".....N/C  
 SPCL-001 \_\_\_\_\_" Custom Platform Width.....\$240.00

√ **PLATFORM TO FLOOR HEIGHT**

- √ <sup>4</sup>SPFH20 18" Platform to Floor Height (with 20" Wheels).....N/C  
 √ <sup>4</sup>SPFH24 22" Platform to Floor Height (with 24" Wheels).....N/C  
 SPFH-001 \_\_\_\_\_" Custom Platform to Floor Height.....\$370.00

√ **WHEEL OPTIONS**

- <sup>7</sup>SPW20 20" Rear Wheels.....N/C  
 <sup>7</sup>SPR24 24 Rear Wheels.....N/C  
 SP-FF Flat Free Wheel Option.....\$125.00

√ **PLATFORM SIDES (Dimension SH)**

- <sup>5</sup>SPS-001 8" Fixed Platform Sides (Foot & Two Sides).....\$500.00  
 <sup>5</sup>SPSC-001 Fixed Sides Custom Height \_\_\_\_\_" (EACH)  
 # \_\_\_\_\_ of sides x \$250.00.....\$ \_\_\_\_\_  
 <sup>6</sup>SPS-002 8" Flip Down Hinged Sides (EACH)  
 # \_\_\_\_\_ of sides x \$300.00.....\$ \_\_\_\_\_  
 <sup>6</sup>SPSC-002 Flip Down Hinged Sides Custom height \_\_\_\_\_" (EACH)  
 # \_\_\_\_\_ of sides x \$350.00.....\$ \_\_\_\_\_  
 <sup>6</sup>SPS-003 8" Adjustable Angle Positive Locking Sides  
 (EACH) # \_\_\_\_\_ of sides x \$400.00.....\$ \_\_\_\_\_  
 <sup>6</sup>SPSC-003 Adjustable Angle Positive Locking Sides Custom height \_\_\_\_\_"  
 (EACH) # \_\_\_\_\_ of sides x \$450.00.....\$ \_\_\_\_\_

#### CUSTOM OPTIONS

- SPCO-001 Tapered Corners (Foot End).....\$250.00  
 SPCO-002 IV Pole & Holder.....\$200.00

#### BASE COVER OPTIONS


- SPBC-001 Custom Base Padding.....Quote Only

#### <sup>3</sup>SIDE PADDING OPTIONS

- <sup>5</sup>SPSP-001 1/2" Poly Side Pads  
 (EACH) # \_\_\_\_\_ of sides x \$85.00.....\$ \_\_\_\_\_  
 <sup>5</sup>SPSP-002 1/2" Sunmate Side Pads    Soft    Medium  
 (EACH) # \_\_\_\_\_ of sides x \$100.00.....\$ \_\_\_\_\_  
 <sup>5</sup>SPSP-003 1" Poly Side Pads  
 (EACH) # \_\_\_\_\_ of sides x \$150.00.....\$ \_\_\_\_\_  
 <sup>5</sup>SPSP-004 1" Sunmate Side Pads    Soft    Medium  
 (EACH) # \_\_\_\_\_ of sides x \$200.00.....\$ \_\_\_\_\_  
 <sup>5</sup>SPSP-005 Vinyl Padding Cover with Zipper.....N/C  
 <sup>5</sup>SPSP-006 Darlex Padding Cover with Zipper  
 (EACH) # \_\_\_\_\_ of sides x \$300.00.....\$ \_\_\_\_\_  
 <sup>5</sup>SPSP-007 Rubatex Padding Cover with Zipper  
 (EACH) # \_\_\_\_\_ of sides x \$200.00.....\$ \_\_\_\_\_  
 <sup>5</sup>SPSP-008 Custom Side Padding.....Quote Only

#### POWER OPTIONS

- SPPO-001 E-Fix Power Drive.....\$7500.00  
 SPPO-002 Battery Box or Vent Tray.....\$299.00  
 Length \_\_\_\_\_" Width \_\_\_\_\_" Height \_\_\_\_\_"  
 SPG-OAD20 One Arm Drive 20" Wheel  
 (choose)  Left    Right.....\$1,500.00  
 SPG-OAD24 One Arm Drive 24" Wheel  
 (choose)  Left    Right.....\$2,500.00



\*\*\*SIDES ARE OPTIONAL\*\*\*

FOOTNOTES ON PAGE 2

**Rehab Series Self Propelling Gurney Option****FOOTNOTES**

*When filling out the measurement form keep in mind certain dimensions.*

## 1: Inside Width: (IDW)

- Without Sides Inside Width = Outside Width of Top Platform
- With Sides With Padding the Width will be Inside the Padding
- With Sides Without Padding the Width will be Inside the ABS

## 1: Inside Length: (IDL)

- Without Sides Inside Length = Outside Length of Top Platform
- With Sides With Padding the Length will be Inside the Padding
- With Sides Without Padding the Length will be Inside the ABS

*Below are equations for figuring out the maximum allowable overall width*

## 2: O.D. Width: (ODW)

- Without Sides same as Inside Width
- With Fixed Sides With Padding:  $1.75'' + (\text{Padding Thickness} \times 2) + \text{Inside Width}$
- With Fixed Sides Without Padding:  $1.75'' + \text{Inside Width}$
- With Hinged Flip Down Sides With Padding:  $2.75'' + (\text{Padding Thickness} \times 2) + \text{Inside Width}$
- With Hinged Flip Down Sides Without Padding:  $2.75'' + \text{Inside Width}$
- With Angle Adjustable Locking Sides With Padding:  $3.75'' + (\text{Padding Thickness} \times 2) + \text{Inside Width}$
- With Angle Adjustable Locking Sides Without Padding:  $3.75'' + \text{Inside Width}$

## 3: Padding Thickness: Thickness of Both Sides of the Platform

## 4: Platform to Floor Height: Measured from the Ground to the Top of the Platform below the ABS

## 5: Maximum Fixed Sides is 3: foot, left side, and/or right side

## 6: Maximum Flip Down or Angle Adjustable Sides 2: left side and/or right side

## 7: Outside Platform Width between the Wheels on the Base: 30'' with Handrims.

## Finished Product Dimensions:

IDW=\_\_\_\_\_inches Inside Width

IDL=\_\_\_\_\_inches Inside Length

ODW=\_\_\_\_\_inches Outside Width

ODL=\_\_\_\_\_inches Outside Length

SH=\_\_\_\_\_inches Side Height

H=\_\_\_\_\_inches Distance from top of steel to floor\*

W=\_\_\_\_\_inches Maximum allowable cart width – for going thru doors

Fits through normal doorways Yes No \*

NOTE: If no pads are ordered, finish product dimensions are inside dimensions of side walls.

Base pad thickness - \_\_\_\_\_inches.

Side pad thickness-\_\_\_\_\_inches.

Hinged sides: right\_\_\_\_\_, left\_\_\_\_\_ or both\_\_\_\_\_

The above dimensions are correct:

\_\_\_\_\_Date:\_\_\_\_\_

Signed

\_\_\_\_\_  
PrintName