

Kidster

Seat to Floor Height Ranges

	Casters Rear Wheels	3 x 1 (K-3X1 OR K- 3X1L)	4 x 1 (K-4X1)	5 1/2 x 1 1/2 (K-515X115)	6 x 2 (K-6X2PN OR K- 6X2SPN)	8 x 2 (K-8X2PN OR K- 8X2SPN)
Rear Wheel Drive	12" Wheels (K-12)	13" - 19"	13" - 19"	13" - 19"	15" - 19"	15" - 19"
	20" Wheels (K-20)	13" - 19"	13" - 19"	13" - 19"	15" - 19"	15" - 19"
	24" Wheels (K-24)	15" - 19"	15" - 19"	15" - 19"	15" - 19"	15" - 19"
Reverse Wheel Drive	12" Wheels (K-12)	13" - 19"	13" - 19"	13" - 19"	15" - 19"	15" - 19"
	20" Wheels (K-20)	15" - 19"	15" - 19"	15" - 19"	15" - 19"	15" - 19"
	24" Wheels (K-24)	15" - 19"	15" - 19"	15" - 19"	15" - 19"	15" - 19"

Special Note: Widened Wheel Base with Armrests off Uprights on 24" Rear Tires

13" & 14" Heights Must Choose Armrests off Uprights

Custom Heights are Available Please Ask!

Rehab TNT

Seat to Floor Height Ranges

	Casters Rear Wheels	5 1/2 x 1 1/2 (K- 515X115)	6 x 2 (K-6X2PN OR K- 6X2SPN)	8 x 2 (K-8X2PN OR K- 8X2SPN)
Rear Wheel Drive	12" Wheels (RT-12)	13" - 19"	15" - 19"	15" - 19"
	20" Wheels (RT-20)	13" - 19"	15" - 19"	15" - 19"
	24" Wheels (RT-24)	15" - 19"	15" - 19"	15" - 19"
Reverse Wheel Drive	12" Wheels (RT-12)	13" - 19"	15" - 19"	15" - 19"
	20" Wheels (RT-20)	15" - 19"	15" - 19"	15" - 19"
	24" Wheels (RT-24)	15" - 19"	15" - 19"	15" - 19"

Special Note: Widened Wheel Base with Armrests off Uprights on 24" Rear Tires

13" & 14" Heights Must Choose Armrests off Uprights

Custom Heights are Available Please Ask!